



# **PRODUCT** SPECIFICATION ENVIRONMENTAL LABEL STERILE

# **PRODUCT INFORMATION**

### PRODUCT DESCRIPTION

Sterile self-adhesive labels (20 labels per sheet, 10 sheets per pouch, 10 pouches per carton (2000 labels per carton total))

#### **PRODUCT FEATURES** Materials

White uncoated paper 70 gsm Acrylic based adhesive Individual label dimensions 45mm x 12mm Sheet size 150mm x 105mm

# Medical Device Classification

Not a medical device (ancillary product)

# MANUFACTURING

30

Celfix, UK

#### PACKAGING

Four layer laminate Easy tear notch Alcohol and VHP compatible Validated three year shelf-life

#### **STERILISATION DETAILS**

Gamma irradiation min. 25kGy Complies with applicable parts of ISO 11137 series

# QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

**SHELF LIFE** 3 years from date of manufacture

#### STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources



# **TECHNICAL INFORMATION**

Compatible with gamma irradiation BSE/TSE free Phthalate free Latex free



# **ORDERING/PACKING INFORMATION**

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ENV200CC	Environmental Label Sterile	n/a	10	200	290 x 170 x 200mm	0.60 kg
ENV200TB	Environmental Label Sterile, Triple-Wrapped	n/a	10	200	290 x 170 x 200mm	0.60 kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.

