



# PRODUCT SPECIFICATION

## ENVIRONMENTAL LABEL

### STERILE

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Sterile self-adhesive labels (20 labels per sheet, 10 sheets per pouch, 10 pouches per carton (2000 labels per carton total))

##### PRODUCT FEATURES

###### Materials

White uncoated paper 70 gsm

Acrylic based adhesive

Individual label dimensions 45mm x 12mm

Sheet size 150mm x 105mm

##### Medical Device Classification

Not a medical device (ancillary product)

##### MANUFACTURING

Celfix, UK

##### PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

Validated three year shelf-life

##### STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of

ISO 11137 series

##### QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

##### SHELF LIFE

3 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
BSE/TSE free  
Phthalate free  
Latex free



## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ENV200CC	Environmental Label Sterile	n/a	10	200	290 x 170 x 200mm	0.60 KG