



# PRODUCT SPECIFICATION

## LARGE (NON-VENTED) DISPENSING PINS PACK STERILE PACKS

### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

Non-vented dispensing pin with safesite valve with luer lock connector

#### PRODUCT FEATURES

##### Materials

Valve: Bayer makrolon polycarbonate

Spike: Bayer lustran ABS

##### Connection type

Luer lock

##### Medical device classification

Not a medical device (pharmacy compounder)

#### MANUFACTURING

B Braun Medical Inc, USA

#### PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

Validated three year shelf-life

#### STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of

ISO 11137 series

#### QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

#### SHELF LIFE

3 years from date of manufacture

#### STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
BSE/TSE free  
Phthalate free  
Latex free



## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
LVDP10CC	B Braun Large Volume Dispensing Pin Sterile 10 Packs of 10	n/a	10	100	290 x 170 x 400mm	1,20 KG