



# PRODUCT SPECIFICATION

## NON-VENTING DISPENSING PINS PACKS STERILE PACKS

### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

Non-vented dispensing pin with safesite valve with luer-lock connector.

#### PRODUCT FEATURES

##### Materials

Valve: Bayer makrolon polycarbonate.  
Spike: Bayer lustran ABS.

##### Connection type

Luer lock.

##### Medical device classification

Not a medical device (pharmacy compounder).

#### MANUFACTURING

B Braun Medical Inc, USA.

#### PACKAGING

Four layer laminate.  
Easy tear notch.  
Alcohol and VHP compatible.  
Validated three year shelf-life.

#### STERILISATION DETAILS

Gamma irradiation min. 25kGy.  
Complies with applicable parts of ISO 11137 series.

#### QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system.

#### SHELF LIFE

3 years from date of manufacture.

#### STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources.



## TECHNICAL INFORMATION

Compatible with gamma irradiation.  
BSE/TSE - considered not to present any risk.\*  
Phthalate free.  
Latex free.

\*This product contains materials of animal origin that are limited to tallow or tallow derivatives. The materials are rigorously processed according to prevailing guidance and standards in order to reduce or eliminate the risk of transmission of BSE/TSE.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
LVDP10CC	Large volume dispensing pins, 10 pack	n/a	10	100	290 x 170 x 400mm	1.20 KG