



PRODUCT SPECIFICATION

SELF-ADHESIVE LABELS STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Sterile self-adhesive labels

PRODUCT FEATURES

Materials

White uncoated paper 70 gsm
Acrylic based adhesive

Medical Device Classification

Not a medical device (ancillary product)

MANUFACTURING

Celfix, UK

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible
Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of
ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016
quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free



ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ENV200CC	Self-adhesive labels 20 labels per sheet	n/a	10	200	290 x 170 x 200mm	0.60 KG