



PRODUCT SPECIFICATION

DELTA PLUS SESSIONAL POLY-CELLULOSE C-FOLDED POUCH WIPES STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Delta Plus is a sterile, low-foaming, neutral detergent formulated for use in cleanrooms and other controlled environments.

Manufactured exclusively from naturally derived ingredients and blended with Water for Injection (WFI), it contains no petrochemical components and carries no hazard pictograms, signal words or hazard statements.

Designed to meet the requirements of EU GMP Annex 1 Section 4.33, Delta Plus provides effective removal of surface contamination, disinfectant residues and bioburden, offering a sustainable and user-friendly solution for the first stage of a validated two-step cleaning and disinfection programmes.

PRODUCT FEATURES

Non-woven polycellulose wipe with sealed edges.

Pre-saturated with Delta Plus to ensure optimum cleaning.

Supplied in a small pack size to minimise waste.

DIMENSIONS

Wipe Size: 23cm x 23cm (9" x 9")

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and sealed in a Grade B (in use) cleanroom. Packaging is performed in a Grade C (in use) cleanroom.

The disinfectant is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Supplied double bagged.

Secondary Packaging: Cardboard carton.

STERILISATION DETAILS

Sterilised by Gamma Irradiation (minimum 25kGy).

SHELF LIFE

Two years from date of manufacture.

STORAGE CONDITIONS

Store in a cool dry place (< 40°C) away from direct sunlight and fluorescent light.

RECYCLING INFORMATION

Component Name:	Component Material:	Disposal:	Guidance:
Wipes	55% Cellulose and 45% Polyester Blend	Non-recyclable	Dispose Via A Licensed Waste Disposal Contractor
Wipes Pouch	Polyester and Polyethylene blend	Non-recyclable	Dispose Via A Licensed Waste Disposal Contractor
Closure Label	White PP Face Material with a Radiation Cured Peel and UV acrylic adhesive	Non-recyclable	Dispose Via A Licensed Waste Disposal Contractor
Polybag	LDPE	Recyclable	Check With Local Authority, Remove Any Labels
Box	Cardboard	Recyclable	Remove Irradiation Labels
Product Label	Gloss Clear Laminated Polypropylene Label with an Acrylic Water-Borne Adhesive	Non-recyclable	Remove From Container
Box Label	Paper Based with Acrylic Water-Borne Adhesive	Non-recyclable Individually	Leave On Cardboard Box
Batch Label (Vellum Label)	Uncoated Paper with an Acrylic Adhesive	Non-recyclable Individually	Leave On Cardboard Box
Cardboard Reinforced Tape	Cardboard	Recyclable	Leave On Cardboard Box

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Pack Qty	Case Qty
	WSNPC1510	Delta Plus Sessional Poly-cellulose C-folded Pouch Wipes Sterile	230 x 230mm	15 wipes per Pouch	15	10 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.