



PRODUCT SPECIFICATION

WHITE SYRINGE CAPS MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Luer lock syringe caps

PRODUCT FEATURES

Materials
DOW HDPE 25055E

Connection type

Luer-Lock (central)

Medical device classification

Class I (sterile)

MANUFACTURING

Baxter Healthcare Ltd UK

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
CAPS10CC	Syringe Cap Sterile 20 Packs of 10	n/a	10	100	280 x 170 x 200mm	0.60 KG
CAPS50CC	Syringe Cap Sterile 10 Packs of 50	n/a	50	500	280 x 170 x 200mm	0.50 KG